

## 检验报告 Certificate of Analysis

No. COA20241214156

免疫球蛋白M(IgM)测定试剂盒(免疫透射比浊法

产品名称

**Product Name** Immunoglobulin M Kit (Turbidimetry

Method)

105-000843-00

批号 144224009

Lot No.

物料编码 **Material Code** 

生产日期 Manufacture

Date

N/A

型号

REF

**Model of the Reagent** 

包装规格

**Package Size** 

有效期

 $R1:1\times 40 \ mL+R2:1\times 10 \ mL$ 

N/A

N/A

2026-09-23 **Expiry date** 

标准依据 产品检验规范 Specification for Product Inspection Standard basis

检测项目	接受标准	检测结果	结果判定
<b>Subjects</b> 1、外观性状 Appearance	Accept criteria	Results	Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	符合/Pass
2、装量 Volume			
	$R2\!\geqslant\!10\text{mL}$	11.03 mL	
净含量 Volume			符合/Pass
	R1≥40mL	41.04 mL	

## 3、性能 ${\tt Performance}$

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within  $\pm 10.0\%$ .

Calibrator: -0.2%

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 0.735(0.570~0.

准确度 900) g/L Accuracy

Testing results of Quality control should be within target range.

测定质控血清,测定结果在靶值范围内。 靶值范围/target range: 1.06(0.82~1.30)

Testing results of Quality control should be within target range.

Control(level 1) :0.78 g/L

符合/Pass

Control(level 2):1.06 g/L

分析灵敏度

吸光度变化△A ≥0.42A (2g/L)。 Analysis of

Sensitivity

Change of absorbance  $\triangle \text{A}\!\geqslant\!0.42\text{A}$  (2g/L)  $^{0.43}$  A

符合/Pass

空白吸光度

试剂空白吸光度<0.3A。

Reagent blank absorbance should be < 0. 0.1 A Blank

符合/Pass

Absorbance

试剂盒在(0.05~4.8)g/L范围内,其回归系数

r ≥0.9900。

样本浓度≤2g/L时,线性绝对偏差≤±0.2g/L; 样本浓度>2g/L时,线性相对偏差≤±10.0%。

线性范围

In the range of  $0.05^4$ . 8g/L, the Linearity Range linear correlation coefficient  $r\!\geqslant\!0.$ 

The absolute deviation should be  $\leq \pm 0.2$ g/L(analyte concentrate  $\leq 2g/L$ ; the relative deviation should be  $\leq \pm 10.0\%$  (

analyte concentrate > 2 g/L)

相关系数R:0.9966;绝对偏差/

Absolute deviation: 0.0;相对 偏差/Relative deviation: Ma

x:6.22%, Min:-4.37%

重复性

批内变异系数 CV≤ 4.0%。

Precision Within-run imprecision CV  $\leq$  4.0%

Within-run

SD:0.0;CV:1.1%

符合/Pass

符合/Pass

结论

Summary

合格/PASS

备注

Remarks

检验: 徐江明Xu Analysis Jiangming | Xu Jiang ming 批准 ( 检验工程师 ): 黄小凤Hua Approved by(Quality Engineer) Xiaofeng

黄小凤Huang

to have thong diastory

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





日期:
Date
日期:
2024-12-14
Date
2024-12-14

## 深圳迈瑞生物医疗电子股份有限公司

深州記論土物医プモナ版の日報公司 SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大夏 邮編: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

