

## 检验报告 Certificate of Analysis

No. COA20250222093

D-二聚体(D-Dimer)测定试剂盒(胶乳增强免疫

产品名称 透射比浊法)

**Product Name** D-Dimer Kit (Particle-enhanced

Immunoturbidimetric Assay Method)

批号 퓆믁 045725003

Lot No.

物料编码 105-004498-00 **Material Code** 

生产日期

Standard basis

Manufacture N/A

Date 标准依据

**REF** 

N/A

N/A

**Model of the Reagent** 

包装规格 **Package Size**  R1:  $1 \times 40$  mL R2:  $1 \times 15$  mL Calibrator:  $6 \times 0.5$  mL

有效期

2026-05-17 **Expiry date** 

产品检验规范 Specification for Product Inspection

检测项目	接受标准	检测结果	结果判定
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance 2、装量 Volume	校准品: 冻干粉,复溶后为清澈透明液体,无沉淀、悬浮物和絮状物。 Calibrator: dry powder, clear liquid without deposit, flocculate or suspended substance after redissolved  R1: 澄清透明的液体,无沉淀、悬浮物和絮状物; R2: 白色混浊液体。 R1, Clear liquid without deposit, flocculate or suspended substance; R2, white turbid liquid	Calibrator:校准品: 冻干粉,复溶后为清澈透明液体,无沉淀、悬浮物和絮状物。 Calibrator: dry powder, clear liquid without deposi t, flocculate or suspended substance after redissolved R1/R2:R1: 澄清透明的液体,无沉淀、悬浮物和絮状物; R2: 白色混浊液体。 R1, Clear liquid without deposit, flocculate or suspended substance; R2, white turbid liquid	符合/Pass
净含量 Volume	R2≥15mL	16.14 mL	符合/Pass

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





 $R1 \ge 40 mL$ 

41.16 mL

## 3、性能 Performance

测定质控血清,测定结果在靶值范围内。

靶值范围/target range: 1.52(1.07~1.97)

Testing results of Quality control

Control(level 1) :1.39  $\mu g/$ 

准确度 Accuracy should be within target range. 测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 15.5(13.1~17.9)

 $\mu g/mL$ 

Testing results of Quality control should be within target range.

Control (level 2):14.75 µg/

mI.

校准品均一性 Homogeneity of Calibrator

校准品的均一性 CV≤10.0%或SD≤0.07μg/

mL.

Homogeneity of Calibrator CV≤10.0%

or SD $\leq$ 0.07  $\mu$  g/mL.

SD: 0. 04 CV: 0. 41%

符合/Pass

符合/Pass

校准品的HBsAg、HIV抗体、HCV抗体、TP抗体检

测应为阴性。

生物安全性 Biosafety

The test of hbsag, hiv antibody, hcv

antibody and tp antibody of the calibration product should be negative

阴性/negative

符合/Pass

试剂空白吸光度<1.5A。 空白吸光度

Blank Reagent blank absorbance should be < 1. 1.48 A

Absorbance

符合/Pass

试剂盒在 (0.5~48) μg/mL范围内, 其回归系

数r ≥0.9900。

样本浓度≤1.4 µg/mL时,线性绝对偏差≤±0. 14 μg/mL; 样本浓度>1.4μg/mL时,线性相

对偏差≤±10.0%。

Linearity Range 9900.

In the range of 0.5~48  $\,\mu\,g/mL$  the linear correlation coefficient  $r \ge 0$ .

相关系数R:0.9980;相对偏差/ Relative deviation: Max:3.4 4%, Min:-3.55%

符合/Pass

The absolute deviation should be  $\leq \pm 0$ . 14 μg/mL(analyte concentrate ≤1.4 μ g/mL); the relative deviation should be  $\leq \pm 10.0\%$  (analyte concentrate >1.4

 $\mu \, g/mL)$ 

重复性

线性范围

批内变异系数 CV≤ 10.0%。

Precision Within-run

Within-run imprecision CV≤10.0%

SD: 0.1; CV: 0.6%

符合/Pass

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





to A Huang Lianforg

**结论**Summary

合格/PASS

备注 Remarks

.....

**检验:** 徐江明Xu **Analysis** Jiangming

**批准 ( 检验工程师 ):** 黄小凤Huang **Approved by(Quality Engineer)** Xiaofeng

2025-02-22

日期:
Date

日期:
Date

was called the control of the called the ca