

检验报告 Certificate of Analysis

No. COA20250306056

产品名称 **Product Name**

糖化血红蛋白质控品 HbA1c Control N

N/A

批号

052425003

Lot No.

N/A

物料编码 **Material Code**

105-002140-00

Package Size

 4×1 mL

生产日期

Manufacture

有效期

REF

型号

N/A

Expiry date

2026-08-26

Date

标准依据 Standard basis

产品检验规范 Specification for Product Inspection

Model of the Reagent

检测项目	接受标准		结果判定
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
包装检验 Package	分装瓶为无色透明硬质玻璃瓶,盖有橡胶瓶内塞,外裹铝制的瓶颈固定;盒贴、瓶贴、标签标识、说明书应完整、清晰。橡胶瓶内塞与玻璃配合后应无泄漏;外裹铝制的瓶颈应无明显划痕、崩缺。 The bottle should be a transparent hard glass bottle with a rubber inner plug and covered with an aluminum bottle neck outside; the packet sticker, bottle sticker, label identification and instructions should be complete and clear. Make sure there is no leakage after the rubber inner plug fits the glass bottle, and there is no obvious scratches or damage to the outer aluminium bottle neck.	分装瓶为无色透明硬质玻璃瓶,盖有橡胶瓶内塞,外裹铝制的瓶颈固定;盒贴、瓶贴、标签标识、说明书完整、清晰。橡胶瓶内塞与玻璃配合后无泄漏;外裹铝制的瓶颈无明显划痕、崩缺。 The bottle were a transparent hard glass bottle with a rubber inner plug and covered with an aluminum bottle neck outside; the packet sticker, bottle sticker, label identification and instructions were complete and clear. Make sure there is no leakage after the rubber inner plug fits the glass bottle, and there is no obvious scratches or damage to the outer aluminium bottle neck.	符合/Pass
参考值表 Reference value table	质控品应提供本批产品的参考值和参考范围。 Provide reference values and reference ranges for this batch of products.	质控品提供本批产品的参考值和参 考范围。 Provide reference values and reference ranges for this batch of products.	符合/Pass

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





瓶内为褐红色冻干粉末状,复溶后应为褐红色液

体, 无沉淀、无悬浮物和絮状物。

颜色性状 Apparent properties Original:Reddish-brown powder; Reconstituted:Reddish-brown liquid without precipitate, suspension and

为褐红色液体, 无沉淀、无悬浮物 和絮状物。 Original:Reddishbrown powder; Reconstituted: Reddish-brown liquid without precipitate, suspension and floccule

瓶内为褐红色冻干粉末状, 复溶后

符合/Pass

2、装量 Volume

每瓶质控品质量(含瓶重)均应不少于5.5g

The weight of each bottle of control (装量 Volume including boottle weight) should be

not less than 5.5g.

6.1 g 符合/Pass

3、性能 Performance

> Hb: 质控品测定结果应在给定的参考范围内。 靶值范围/target range: 126(114~138)

umo1/L

floccule

 $Hb:125.43 \mu mol/L$

Testing results of Quality control should be within target range.

HbA1c%: 质控品测定结果应在给定的参考范围 参考值及参考范围 内。 靶值范围/target range: 5.43%(4.

Reference and

89% 5.97%) reference range Testing results of Quality control

should be within target range.

HbA1c: 质控品测定结果应在给定的参考范围 内。 靶值范围/target range: 4.50(4.05 ~4.95) µ mo1/L

Testing results of Quality control should be within target range.

HbA1c%: 5.42%

 $HbA1c:4.48 \mu mo1/L$

Hb: 瓶内SD/ SD within<6.40 μmol/L, 瓶

内CV%/ CV% within<5.00%

Hb:SD:0.37;CV:0.30%

瓶内均一性 Inter-vial homogeneity HbA1c%: 瓶内SD/ SD within<0.50 %, 瓶内

CV%/ CV% within<5.00%

HbA1c%: SD: 0.04%; CV: 0.69%

符合/Pass

符合/Pass

HbA1c: 瓶内SD/ SD within<0.55 μmol/L或

瓶内CV%/ CV% within<5.00%

HbA1c:SD:0.04;CV:0.96%

HBsAg、HIV抗体、HCV抗体、梅毒螺旋体TP抗体

检测应为阴性。

生物安全性 Biosafety

HBsAg , HIV antibody, HCV antibody, and 阴性/negative

TP antibody testing results of the

control should be negative.

符合/Pass

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





to A Huang Lianforg

结论Summary

合格/PASS

备注

Remarks

检验: 徐江明Xu **Analysis** Jiangming | Xu Jiang ming 批准 (检验工程师): 黄小凤Huang Approved by(Quality Engineer) Xiaofeng

日期:
Date
日期:
2025-03-06
Date
2025-03-06

深圳迈瑞生物医疗电子股份有限公司
SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680
Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com