

检验报告 Certificate of Analysis

No. COA20250325026

产品名称 **Product Name**

肌酸激酶同工酶MB校准品 CK-MB Calibrators

N/A

批号

Lot No.

2025020101

N/A Model of the Reagent

物料编码

105-005913-00 **Material Code**

 $C0:1\times 2.0$ mL, $C1:1\times 2.0$ mL, $C2:1\times 2.$

Package Size

生产日期

2025-03-24

有效期

REF

型号

2026-09-23 **Expiry date**

Date

标准依据 Standard basis

Manufacture

产品检验规范 Specification for Product Inspection

检测项目 检测结果 结果判定 接受标准 Subjects Accept criteria Results Conclusion

1、外观性状 Appearance

> 校准品复融前外观呈现粉末块状固体,不起泡, 不塌陷,复融后为清澈透明液体,无沉淀、无悬 浮物、无絮状物; 分装瓶为棕色的玻璃瓶, 盖有 塑料外盖; 盒贴、瓶贴、标签标识、说明书、校 准卡完整、清晰牢固;标贴信息正确、完整;塑 料外盖与玻璃瓶身配合后无泄漏; 塑料外盖无明 显划痕、崩缺。

The calibrators are solid power without bubble or collapse before dissolved. And are clear liquid free of precipitate, suspension and

外观和性状 Appearance floccules after dissolved.

符合标准/Qualified

符合/Pass

The containers are brownish glass bottles with tightly closed caps; The whole package, including labels for boxes and bottles, pack inserts and calibration cards, is complete and clear; The content of labels is correct and complete. The containers are without leakage and obvious defect

or breakdown.

3、性能 Performance

瓶内均一性 Uniformity in bottle

 $SD \leq 0.28 ng/mL$ $SD \leq 0.28 ng/mL$

C0:0.00 ng/mL

符合/Pass

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





CV ≤ 8.0% CV ≤ 8.0%

CV≤8.0% CV≤8.0%

偏倚在±10.0%范围内。

The Bias of Trueness Control Sample is H:3.0%

less than $\pm 10.0\%$.

赋值准确性 符合/Pass

偏倚在土10.0%范围内。

The Bias of Trueness Control Sample is $\ L{:}1.4\%$

less than $\pm 10.0\%$.

结论 Summary 合格/PASS

备注 Remarks

检验: 邓月Deng Yue スプリープログイル 批准(检验工程师): Approved by(Quality Engineer) 林朋Lin Peng

日期:
Date
日期:
2025-03-25
Date
2025-03-25

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680
Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com



林朋 Lingeng