

## 检验报告 Certificate of Analysis

No. COA20250109032

产品名称 **Product Name**  类风湿因子(RF)测定试剂盒(免疫比浊法)

Rheumatoid Factor Kit (

Immunoturbidimetric Method)

批号

147025001 Lot No.

物料编码

**Material Code** 

生产日期

Manufacture

Date 标准依据

N/A

105-007676-00

REF

N/A

N/A

型号

Model of the Reagent

包装规格

Package Size

 $R1:1\times40$  mL+R2:1×11 mL+Calibrato

 $r:5\times0.5$  mL

有效期

**Expiry date** 

2026-02-20

产品检验规范 Specification for Product Inspection Standard basis

检测项目	接受标准		结果判定
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	渗漏,标签完好清晰,是无标签脱	符合/Pass
外观物理性状 Appearance	校准品: 无色或浅黄色液体。 Calibrator: Colorless or light yellow liquid.	Calibrator:校准品: 无色或浅黄 色液体。 Calibrator: Colorless or light yellow liquid.	符合/Pass
	澄清透明的液体,无沉淀、悬浮物和絮状物。	R1/R2:澄清透明的液体,无沉淀、 悬浮物和絮状物。	
	Reagent: Clear liquid without deposition and floccule	Reagent: Clear liquid without deposition and floccule	
2、装量 Volume			
	R2≥11mL	12.02 mL	
净含量 Volume			符合/Pass
	R1≥40mL	41.10 mL	

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





Calibrator(5 level)≥0.5mL

a: 0.50mL;b: 0.50mL;c: 0.50mL;d: 0.50mL;e: 0.50mL

## 3、性能 Performance

Accuracy

测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within  $\pm 10.0\%$ .

b:7.1%

测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within  $\pm 10.0\%$ .

c:1.2%

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 31.3(20.2~42.4)

TU/mL

Control(level 1) :30.6 IU/mL

Testing results of Quality control 准确度 should be within target range.

测定质控血清,测定结果在靶值范围内。

靶值范围/target range: 72.5(47.0~98.0)

IU/mL

Control(level 2):69.2 IU/mL

Testing results of Quality control should be within target range.

测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should  $% \left( 1\right) =\left( 1\right) \left( 1\right) \left$ 

be within  $\pm 10.0\%$ .

d:0.4%

测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should  $% \left\{ 1\right\} =\left\{ 1$ 

be within  $\pm 10.0\%$ .

e:0.7%

分析灵敏度 吸光度变化△A ≥0.01A ( 20 IU/mL)。

Analysis of Change

Change of absorbance  $\triangle A \geqslant 0.01A$  (20 I 0.03 A

U/mL)

符合/Pass

符合/Pass

校准品均一性

Sensitivity

Homogeneity of

Calibrator

校准品的均一性 CV≤8.0%。

Homogeneity of Calibrator CV≤8.0%.

1.6%

符合/Pass

校准品的HBsAg、HIV抗体、HCV抗体、TP抗体检

测应为阴性。

生物安全性 The test of hbsag, hiv antibody, hcv

Biosafety antibody and tp antibody of the

calibration product should be negative

阴性/negative

符合/Pass

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680
Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





试剂空白吸光度<1.0A。 空白吸光度

Reagent blank absorbance should be<1. 0.1 A Blank

符合/Pass

Absorbance

试剂盒在(10~500) IU/mL范围内, 其回归系数

r ≥0.9900。

样本浓度≤40 IU/mL时,线性绝对偏差≤±4 IU/mL; 样本浓度>40 IU/mL时,线性相对偏差

线性范围

In the range of  $10^500 \text{ IU/mL}$ , the

linear correlation coefficient  $r \ge 0$ . Linearity Range 9900.

> The absolute deviation should be  $\leq$   $\pm 4$ IU/mL(analyte concentrate ≤40 IU/m

L); the relative deviation should be  $\leq$  $\pm 10.0\%$  (analyte concentrate >40 IU/mL

重复性

批内变异系数 CV≤ 8.0%。

Precision Within-run imprecision CV≤8.0% Within-run

SD: 0.7; CV: 2.4%

0%, Min:-4.68%

相关系数R:0.9981;相对偏差/

Relative deviation: Max:4.5

符合/Pass

to Not Huang Lianforg

符合/Pass

结论

Summary

合格/PASS

备注

Remarks

徐江明Xu Analysis Jiangming

日期:

2025-01-08 Date

| Xu Jiang ming | 批准(检验工程师): 黄小凤Hua | Approved by(Quality Engineer) Xiaofeng

日期:

Date

2025-01-09

黄小凤Huang



SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

