

N/A

N/A

## 检验报告 Certificate of Analysis

No. COA20241009054

免疫球蛋白A(IgA)测定试剂盒(免疫透射比浊法

产品名称

**Product Name** Immunoglobulin A Kit (Turbidimetry

Method)

批号

型号 148824005 Lot No. **Model of the Reagent** 

物料编码 包装规格 105-004618-00

 $R1:1\times 36 \text{ mL+R2:}1\times 18 \text{ mL}$ **Material Code Package Size** 

REF

生产日期

有效期 Manufacture 2026-08-20 N/A**Expiry date** 

Date 标准依据

产品检验规范 Specification for Product Inspection Standard basis

检测项目	接受标准	检测结果	
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	符合/Pass
2、装量 Volume			
净含量 Volume	R2≥18mL	18.77 mL	符合/Pass
	R1≥36mL	36.80 mL	

## 3、性能 Performance

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within  $\pm 10.0\%$ .

Calibrator: 0.6%

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 1.67( 1.28~2.06)

准确度 Accuracy

Control(level 1) :1.68 g/L

符合/Pass

Testing results of Quality control should be within target range.

测定质控血清,测定结果在靶值范围内。

靶值范围/target range: 2.47(1.90~3.04)

g/L

Control(level 2):2.52 g/L

Testing results of Quality control should be within target range.

分析灵敏度

吸光度变化△A ≥0.63A (2.5g/L)。

Change of absorbance  $\triangle A \ge 0.63A$  (2.5g/ 0.86 A Analysis of

Sensitivity

符合/Pass

空白吸光度

试剂空白吸光度<0.3A。

Reagent blank absorbance should be < 0. 0.1 A Blank Absorbance

符合/Pass

试剂盒在 (0.2~8.0) g/L范围内, 其回归系数r

≥0.9900。

样本浓度≤2.4g/L时,线性绝对偏差≤±0.24 g/L;样本浓度>2.4g/L时,线性相对偏差≤±1

0.0%

线性范围

In the range of  $0.2^8.0$  g/L, the

Linearity Range linear correlation coefficient  $r\!\geqslant\!0.$ 

相关系数R:0.9954;绝对偏差/ Absolute deviation: -0.08;相 对偏差/Relative deviation: Max:8.94%, Min:-4.52%

符合/Pass

The absolute deviation should be  $\leq \pm 0$ . 24 g/L (analyte concentrate  $\leq 2.4$ g/L); the relative deviation should be  $\leq \pm 1$ 0.0%(analyte concentrate > 2.4 g/L)

重复性

批内变异系数 CV≤ 4.0%。

Precision Within-run imprecision CV  $\leq$  4.0% Within-run

SD: 0.0; CV: 0.7%

符合/Pass

结论

Summary

合格/PASS

备注

Remarks

检验: 徐江明Xu Analysis Jiangming | Xu Jiang ming 批准 ( 检验工程师 ): 黄小凤Hua Approved by(Quality Engineer) Xiaofeng

黄小凤Huang

to have thong diastory

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





日期: 日期: 2024-10-09 2024-10-09 Date Date

## 深圳迈瑞生物医疗电子股份有限公司

深州記論土物医プモナ版の日報公司 SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大夏 邮編: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

