

检验报告 Certificate of Analysis

No. COA20250118073

产品名称 **Product Name** 碱性磷酸酶(ALP)测定试剂盒(AMP缓冲液法) Alkaline Phosphatase Kit (IFCC

N/A

批号

140325001

Modified Method)

型号

REF

Model of the Reagent

N/A

Lot No.

物料编码 105-000855-00 包装规格

Material Code

Package Size

 $R1:6\times40~mL+R2:2\times32~mL$

生产日期

Manufacture

N/A

有效期

2026-07-13

Date 标准依据 **Expiry date**

产品检验规范 Specification for Product Inspection Standard basis

 检测项目	 接受标准		 结果判定
Subjects 1、外观性状 Appearance	Accept criteria		Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	符合/Pass
2、装量 Volume			
净含量	R2≥32mL	33.04 mL	
Volume	R1≥40mL		符合/Pass
	WI > TOWL	40.67 mL	

3、性能 ${\tt Performance}$

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within $\pm 10.0\%$.

Calibrator: 0.5%

测定质控血清, 测定结果在靶值范围内。

准确度 Accuracy

靶值范围/target range: 100(85~115)U/L Testing results of Quality control

should be within target range.

Control(level 1):98.5 U/L

符合/Pass

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 241(205~277)U/L

Testing results of Quality control should be within target range.

Control(level 2):241.2 U/L

分析灵敏度

吸光度变化率 ≥0.014 A/min (75U/L)。

Change of absorbance △A≥0.014 A/min 0.015 A/min Analysis of

Sensitivity (75U/L)

符合/Pass

空白吸光度

试剂空白吸光度<2.5A。

Reagent blank absorbance should be < 2. 0.4 A Blank

Absorbance

符合/Pass

空白吸光度变化率 空白吸光度变化率<0.007A/min

Absorbance Rate Reagent blank absorbance<0.007A/min

0.003 A/min

符合/Pass

试剂盒在(5~800) U/L范围内, 其回归系数r

样本浓度≤120U/L时,线性绝对偏差≤±12.0 U/L;样本浓度>120U/L时,线性相对偏差≤±1

0.0%。

线性范围

In the range of $5^{\sim}800$ U/L, the linear

Linearity Range correlation coefficient $r \ge 0.9900$. The absolute deviation should be \leq \pm

12.0 U/L(analyte concentrate \leq 120 U/ L); the relative deviation should be \leq $\pm 10.0\%$ (analyte concentrate >120 U/L)

相关系数R:1.0000;相对偏差/

Relative deviation: Max:0.9

8%, Min:-0.49%

符合/Pass

重复性

批内变异系数 CV≤3.0%。


Precision Within-run imprecision CV≤3.0% Within-run

SD:1.7;CV:1.7%

符合/Pass

结论 Summary

合格/PASS



SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





备注 Remarks

日期:

Date

日期:

Date

检验: 徐江明Xu Analysis

Jiangming

2025-01-18

| Xu Jiang ming 批准 (检验工程师): 黄小凤Hua Approved by(Quality Engineer) Xiaofeng

黄小凤Huang

2025-01-18

the A Huang Stanforg

深圳迈瑞生物医疗电子股份有限公司

深州記論土物医プモナ版の日報公司 SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大夏 邮編: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

