

检验报告 Certificate of Analysis

No. COA20241221056

C-反应蛋白(CRP)测定试剂盒(免疫透射比浊法) REF 产品名称

C-Reactive Protein Kit (Turbidimetry **Product Name**

Method)

批号 148924021

Lot No.

物料编码 105-000880-00 **Material Code**

生产日期

Manufacture N/A

Date

N/A

N/A

型号

Model of the Reagent

包装规格

Package Size

 $R1:1\times 40 \text{ mL}+R2:1\times 10 \text{ mL}$

有效期

Expiry date

2026-01-24

标准依据 产品检验规范 Specification for Product Inspection Standard basis

检测项目 Subjects 1、外观性状 Appearance	接受标准 Accept criteria	检测结果 Results	结果判定 Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	符合/Pass
2、装量 Volume			
净含量	R2≥10mL	11.26 mL	
77百里 Volume	R1≥40mL	41.10 mL	符合/Pass

3、性能 ${\tt Performance}$

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within $\pm 10.0\%$.

Calibrator: 0.7%

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 17.5(12.1~22.9)

准确度 mg/L

Accuracy Testing results of Quality control should be within target range.

> 测定质控血清,测定结果在靶值范围内。 靶值范围/target range: 42.8(29.9~55.7)

mg/L

Testing results of Quality control should be within target range.

Control(level 1) :18.2 mg/L

符合/Pass

Control(level 2):44.9 mg/L

吸光度变化△A在(0.05~0.50)A (40 mg/ 分析灵敏度

Analysis of Sensitivity A (40 mg/L)

Change of absorbance $\triangle A \ge (0.05^{\circ}0.50)$

0.11 A

空白吸光度

试剂空白吸光度<0.3A。

Reagent blank absorbance should be < 0. Blank

Absorbance

0.0 A

符合/Pass

符合/Pass

试剂盒在(2~250) mg/L范围内, 其回归系数r

≥0.9900。

样本浓度≤60mg/L时,线性绝对偏差≤±6.0m g/L;样本浓度>60mg/L时,线性相对偏差≤±1

0.0%。

线性范围 Linearity Range

In the range of $2^{\sim}250$ mg/L, the linear correlation coefficient r≥0.9900.

The absolute deviation should be $\leq \pm 6$. 0 mg/L(analyte concentrate ≤ 60 mg/L); the relative deviation should be $\leq \pm 1$ 0.0%(analyte concentrate >60 mg/L)

相关系数R:0.9976;绝对偏差/ Absolute deviation: 2.5;相对 偏差/Relative deviation: Ma

x:3.95%, Min:-5.12%

符合/Pass

重复性

批内变异系数 CV≤ 6.0%。

Precision Within-run imprecision CV≤6.0% Within-run

SD: 0.2; CV: 1.0%

符合/Pass

结论

Summary

合格/PASS

备注

Remarks

检验: Jiangming Analysis

| Xu Jiang ming 批准 (检验工程师): 黄小凤Hua Approved by(Quality Engineer) Xiaofeng

the Alwang Lianforg

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





日期: 日期: 2024-12-21 2024-12-21 Date Date

深圳迈瑞生物医疗电子股份有限公司

深州記論土物医プモナ版の日報公司 SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大夏 邮編: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

