

## 检验报告 Certificate of Analysis

No. COA20250412083

产品名称 **Product Name** 

肌酸激酶同工酶校准品 CK-MB Calibrator

N/A

批号

151425004

N/A

Lot No.

物料编码

Model of the Reagent

105-001132-00 **Material Code** 

Package Size

 $3 \times 1$  mL

生产日期

Manufacture Date

N/A

有效期 **Expiry date** 

**REF** 

型号

2026-09-30

标准依据 Standard basis

产品检验规范 Specification for Product Inspection

检测项目	接受标准		结果判定
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
包装检验 Package	分装瓶为棕色硬质玻璃瓶,盖有橡胶瓶内塞和塑料外盖;盒贴、瓶贴、标签标识、说明书应完整、清晰。橡胶瓶内塞与玻璃配合后应无泄漏;塑料外盖应无明显划痕、崩缺。 The bottle should be a brown hard glass bottle with a rubber inner plug and a plastic outer lid; the packet sticker, bottle sticker, label identification, and instructions should be complete and clear. Make sure there is no leakage after the rubber inner plug fit the bottle, and there is no obvious scratches or damages to the plastic outer lid.	分装瓶为棕色硬质玻璃瓶,盖有橡胶瓶内塞和塑料外盖;盒贴、瓶贴、标签标识、说明书完整、清晰。橡胶瓶内塞与玻璃配合后无泄漏;塑料外盖无明显划痕、崩缺。 The bottle were a brown hard glass bottle with a rubber inner plug and a plastic outer lid; the packet sticker, bottle sticker, label identification, and instructions were complete and clear. Make sure there is no leakage after the rubber inner plug fit the bottle, and there is no obvious scratches or damages to the plastic outer lid.	符合/Pass
参考值表 Reference value table	校准品应可提供所包含检测项目的参考值表。 For the calibrator, provide a reference value table for the testing items included.	校准品提供所包含检测项目的参考 值表。 For the calibrator, provide a reference value table for the testing items included.	符合/Pass
颜色性状 Apparent properties	瓶内为白色冻干粉末。复溶后为清澈透明无色液体,无沉淀、无悬浮物和絮状物Original:White powder;Reconstituted:Clear and colorless liquid without precipitate, suspension and floccule	瓶内为白色冻干粉末。复溶后为清澈透明无色液体,无沉淀、无悬浮物和絮状物 Original:White powder;Reconstituted:Clear and colorless liquid without precipitate, suspension and floccule	符合/Pass

## 深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





符合/Pass

符合/Pass

2、装量 Volume

每瓶肌酸激酶同工酶校准品重量(含瓶重)应不

少于6g。

装量 The weight of each bottle of

Volume calibrator (including boottle weight)

should be not less than 6g.

3、性能 Performance

瓶内SD/ SD within<3.00 U/L或瓶内CV%/ C 瓶内均一性

Inter-vial V% within<2.0% SD: 0.45 CV: 0.32% 符合/Pass

8 g

homogeneity

HBsAg、HIV抗体、HCV抗体、梅毒螺旋体TP抗体

检测应为阴性。

生物安全性 HBsAg , HIV antibody, HCV antibody, and 阴性/negative 符合/Pass Biosafety

TP antibody testing results of the calibrator should be negative.

允许偏倚界限值 (绝对偏倚或相对偏倚): ±

2.40 U/L或±10.0%

Absolute bias or relative bias:  $\pm 2.40$ 

U/L or  $\pm 10.0\%$ 

赋值准确性 Accuracy

允许偏倚界限值 (绝对偏倚或相对偏倚): ±

2.40 U/L或±10.0%

Absolute bias or relative bias:  $\pm 2.40$ 

U/L or  $\pm 10.0\%$ 

Reference seruml:绝对偏倚/

Absolute bias:-0.23;相对偏倚/

Relative bias:-1.43%

Reference serum2:绝对偏倚/

Absolute bias:-0.10;相对偏倚/

Relative bias:-0.24%

结论

合格/PASS **Summary** 

备注

检验:

Remarks

徐江明Xu Analysis

Jiangming

日期:

2025-04-12 Date

| Xu Jiang ming | 批准 ( 检验工程师 ): 黄小凤Hua | Approved by(Quality Engineer) | Xiaofeng

Date

2025-04-12

黄小凤Huang



the Huang Lianforg