

检验报告 Certificate of Analysis

No. COA20250306137

产品名称 Product Name 补体因子C3(C3)测定试剂盒(免疫透射比浊法) Complement C3 Kit (Turbidimetry Method

N/A

Product Nai 批号 omplement of kit (luibluimetry method

--, -

Lot No.

143925002

型号 Model of the Reagent

N/A

物料编码 Material Code

105-004617-00

包装规格 Package Size

 $R1:1\times36$ mL+R2:1×18 mL

生产日期

Manufacture N/A

有效期

REF

Date 标准依据 Expiry date

2027-02-18

小座水塘 Standard basis 产品检验规范 Specification for Product Inspection

检测项目	接受标准	检测结果	结果判定
Subjects 1、外观性状 Appearance	Accept criteria	Results	Conclusion
外包装检验 Outer packaging inspection	试剂外包装完整,标识清晰,密封良好,无破裂现象。 试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	试剂外包装完整,标识清晰,密封良好,无破裂现象。试剂瓶液体无渗漏,标签完好清晰,是无标签脱落、污损现象,批号、效期信息完整。 The reagent package is complete, clearly marked, well sealed, and free from cracking.	符合/Pass
外观物理性状 Appearance	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	清澈透明的液体,无沉淀、悬浮物和絮状物。 Clear liquid without deposition and floccule.	符合/Pass
2、装量 Volume			
	R2≥18mL	19.00 mL	
净含量 Volume			符合/Pass
	R1≥36mL	37.03 mL	

3、性能 Performance

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.
地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680
Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





测定校准品,测定结果与靶值的相对偏差在±1

0.0% 以内。

The bias of calibrator testing should

be within $\pm 10.0\%$.

Calibrator: -4.8%

测定质控血清, 测定结果在靶值范围内。

靶值范围/target range: 1.04(0.83~1.25)

准确度 g/LAccuracy

Control(level 1) :1.05 g/L

符合/Pass

should be within target range. 测定质控血清,测定结果在靶值范围内。

Testing results of Quality control

靶值范围/target range: 1.50(1.2~1.8)g/L Control(level 2):1.55 g/L

Testing results of Quality control should be within target range.

分析灵敏度

吸光度变化△A ≥0.37A (1.5g/L)。

Change of absorbance $\triangle A \ge 0.37A$ (1.5g/ 0.60 A Analysis of

Sensitivity L) 符合/Pass

空白吸光度

试剂空白吸光度<0.10A。

Reagent blank absorbance should be < 0. 0.06 A Blank

Absorbance

符合/Pass

试剂盒在(0.04~3.3)g/L范围内,其回归系数

r ≥0.9900。

样本浓度≤1.44g/L时,线性绝对偏差≤±0. 144g/L; 样本浓度>1.44g/L时, 线性相对偏差

 $\leq \pm 10.0\%$ 。

In the range of $0.04^{\circ}3.3g/L$, the 线性范围

linear correlation coefficient $r \ge 0$.

The absolute deviation should be \leq \pm

0.144 g/L(analyte concentrate ≤1.44g/ L); the relative deviation should be \leq $\pm 10.0\%$ (analyte concentrate > 1.44 g/L

Linearity Range 9900.

)

相关系数R:0.9992;绝对偏差/ Absolute deviation: Max:0. 033, Min:-0.038; 相对偏差/ Relative deviation: Max:2.1

5%, Min:-3.31%

符合/Pass

重复性 批内变异系数 CV≤ 3.0%。 Precision

Within-run imprecision CV \leq 3.0% Within-run

SD: 0. 0: CV: 0. 6%

符合/Pass

结论

Summary

合格/PASS

备注

Remarks

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





检验: 徐江明Xu **Analysis**

Jiangming

批准 (检验工程师): 黄小凤Huang **Approved by(Quality Engineer)** Xiaofeng

日期: 2025-03-06 Date

Date

2025-03-06

深圳迈瑞生物医疗电子股份有限公司

深州記論土物医プモナ版の日報公司 SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大夏 邮編: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

