

N/A

N/A

 2×50 tests

2026-05-19

检验报告 Certificate of Analysis

REF

Expiry date

No. COA20250402034

人类免疫缺陷病毒抗原抗体(HIV)测定试剂盒(

产品名称 化学发光免疫分析法)

Product Name Antigen and Antibodies to Human

Immunodeficiency Virus(CLIA)

批号

2024080122 Lot No. Model of the Reagent

物料编码 包装规格

105-004234-00 **Material Code** Package Size

生产日期

Manufacture 2024-11-20

Date

标准依据 产品检验规范 Specification for Product Inspection Standard basis

..... 检测项目 接受标准 检测结果 结果判定 Subjects **Accept criteria** Results **Conclusion** 1、外观性状 Appearance 试剂盒各组分齐全、完整、液体无渗漏; 包装标签清晰、无破损; Ra组分为棕色含固体微粒的液体; Rb、Rc、Rd组分为清澈透明的液体 Reagent Componets are complete, no 外观和性状 符合标准/Qualified 符合/Pass leakage. Labels are clear and well-Appearance stamped. Ra is composed of brown liquid containing solid particles. Rb, Rc and Rd are composed of clear and transparent liquid. 2、装量 Volume Rb≥3.5mL 3.8 mL $Rc \ge 3.5 mL$ 装量 符合/Pass 3.8 mL Volume Ra≥3.8mL 4.0 mL

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





Rd≥5.3mL 5.5 mL

3、性能 Performance

最低检测限 Limit of Detection

最低检测限≤2.0 IU/mL

The minimum detection is ≤2.0 IU/mL

0.5 IU/mL

符合/Pass

符合/Pass

灵敏度参考品(AE重组型)阳性检出率应≥3/5

HIV-1 CRF_AE sensitivity reference

material: ≥3/5

AE重组型

 $(S-AE1 \sim S-AE5) : 4/5$

灵敏度参考品(B/B'亚型)阳性检出率应≥

HIV-1 subtype B/B' sensitivity reference material: $\geq 3/5$

B/B'亚型

 $(S-B1 \sim S-B5) : 5/5$

灵敏度参考品 Analysis Sensitivity Reference Material

灵敏度参考品(BC重组型)阳性检出率应≥3/5

BC重组型 HIV-1 CRF_BC sensitivity reference

 $(S-BC1 \sim S-BC5) : 5/5$

material: $\geq 3/5$ 10份HIV-1 p24抗原灵敏度国家参考品或经国家

参考品标化的灵敏度参考品L10应判断为阴性, L1、L2、L3和L4(最低检出量≤2.5 IU/mL)

应判断为阳性

HIV-1 P24 antigen sensitivity reference material:L1、L2、L3、L4 (the 2.5IU/m1); L10阴性Negetive lowest detection level is ≤2.5 IU/m

L): Positive, L10 is negtive

HIV-1 p24抗原:L1阳性Positiv e; L2阳性Positive; L3阳性 Positive; L4阳性Positive (<

变异系数CV≤ 10%

The coefficient of variation is ≤10% HIV-1阳性:2%

重复性

Repeatability

变异系数CV≤ 10%

The coefficient of variation is ≤10% HIV-2阳性:2%

符合/Pass

变异系数CV≤ 10%

The coefficient of variation is ≤10% HIVP24阳性:2%

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com





符合/Pass

符合/Pass

变异系数CV≤ 10%

The coefficient of variation is ≤10% HIV阴性:3%

10份HIV-1 p24抗原阳性参考品符合率 (+/+)

为10/10

HIV-1 P24 Antigen Positive Reference

Material: 10/10.

HIV-1 p24抗原 (P1-P10):10/10

HIV-1型M组阳性参考品符合率 (+/+) 为14/14

HIV-1 M Group Positive Reference

(P1~P14):14/14

阳性参考品符合率 Material: 14/14.

Positive

reference rate HIV-1型0组阳性参考品中,阳性符合率应≥1/

HIV-1 O Group Positive Reference

Material: ≥1/3

HIV-1型0

HIV-1型M

(P15∼P17) :3/3

HIV-2型阳性参考品中,阳性符合率应≥2/3

HIV-2 Positive Reference Material: ≥

2/3

HIV-2

(P18~P20) :3/3

20份HIV-1 p24抗原阴性参考品符合率 (-/-)

为20/20。

HIV-1 P24 Antigen Negative Reference

HIV-1 p24抗原(N1~N20):20/20

阴性参考品符合率 Material: 20/20.

Negative

reference rate 13份人类免疫缺陷病毒抗体阴性参考品符合

率 (-/-) 应为13/13

HIV-1/2 Antibody Negative Reference

Material:13/13.

HIV-1/2抗体

 $(N1^{\sim}N13):13/13$

结论 **Summary**

合格/PASS

备注

Remarks

检验: 邓月Deng Yue **Analysis**

AH Yeng Xue

批准(检验工程师):

Approved by(Quality Engineer) 林朋Lin Peng

林朋 lingeng

日期:

2025-04-02 Date

日期: Date

2025-04-02

深圳迈瑞生物医疗电子股份有限公司

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. 地址: 深圳市南山区高新技术产业园区科技南十二路迈瑞大厦 邮编: 518057 电话: 0755-8188 8998 传真: 0755-26582680 Adress: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 Website: www.mindray.com

